

Waveform digitizing for LaBr₃/NaI phoswich detector^{*}

SHE Rui(余睿)^{1,2} JI Jian-Feng(季建峰)^{1,2;1)} FENG Hua(冯骅)^{1,2} LI Hong(李红)^{1,2} HAN Dong(韩冬)^{1,2}

¹ Department of Engineering Physics, Tsinghua University, Beijing 100084, China

² Key Laboratory of Particle & Radiation Imaging (Tsinghua University), Ministry of Education, Beijing 100084, China

Abstract: The detection efficiency of phoswich detector starts to decrease when Compton scattering becomes significant. Events with energy deposit in both scintillators, if not rejected, are not useful for spectral analysis as the full energy of the incident photon cannot be reconstructed with conventional readout. We show that once the system response is carefully calibrated, the full energy of those double deposit events can be reconstructed using a waveform digitizer as the readout. Our experiment suggests that the efficiency of a photopeak at 662 keV can be increased by a factor of 2 using our LaBr₃/NaI phoswich detector.

Key words: X-ray detector, phoswich, LaBr₃, waveform digitizing

PACS: 95.55.Ka **DOI:** 10.1088/1674-1137/38/8/086005

1 Introduction

The phoswich detector, a combination of two scintillation crystals with distinct light decay times, has been widely used in hard X-ray astronomy, in particular where large detection areas are needed. In the past, a combination of NaI(Tl) and CsI(Na), with a decay time of 630 ns and 250 ns respectively, has often been used as a phoswich. For example, the Phoswich Detection System (PDS) onboard BeppoSAX [1], which is sensitive in the energy range of 15 to 300 keV, has resulted in fruitful scientific outcome. Recently, LaBr₃ doped with Ce was found to be a new scintillation material with better energy resolution and much shorter decay time (16 ns) than that of NaI [2]. Several studies suggest that a phoswich detector using a combination of LaBr₃ and NaI has great potential for X-ray astronomy [3–5], especially for time domain astronomy and transient detection.

In phoswiches, the top crystal acts as the sensitive volume, while the bottom crystal is used as an active veto to shield background. The sensitive layer should be much thinner than the veto layer in order for effective background rejection. For X-rays with energy above ~200–300 keV, Compton scattering becomes significant in detectors, leading to a low detection efficiency in the phoswich due to energy deposit in the bottom layer by Compton events. Even if not rejected, events with energy deposit in both layers cannot be used for spectroscopy due to different decay time and light yields. This makes it difficult for phoswich detectors to detect X-rays above 300 keV, and interesting science in this energy range

may be missed, such as the Galactic electron-positron annihilation at 511 keV [6], or gamma-ray bursts with high peak energies [7].

In a previous study [5], we have demonstrated the technical readiness of a large area ($\Phi 101.6$ mm) LaBr₃:Ce/NaI(Tl) phoswich detector and its potential capability for hard X-ray astronomy. This phoswich detector shows good low energy response (below 5 keV), fine energy resolution ($\sim 10\%$ at 60 keV), and excellent spatial unification (less than 1% in pulse amplitude across the sensitive surface). In the pulse width spectrum, the two peaks caused by energy deposit in either LaBr₃:Ce or NaI(Tl) are distinctly separated as per their large decay time ratio (250:16).

Here, we show that if a waveform digitizer is used as the readout, Compton events with energy deposit in both layers can be recognized and the full energy of incident photons can be reconstructed by fitting the pulse waveform. This can improve the high energy response of phoswich detectors.

2 Instrument setup

The phoswich detector consists of cylindrical LaBr₃:Ce and NaI(Tl) scintillators with a diameter of 101.6 mm. The top LaBr₃:Ce layer is 6 mm thick and the bottom NaI(Tl) layer is 40 mm thick. With such a design, the detector can achieve a full energy efficiency of 20% up to 300 keV [5]. The entrance window is 0.22 mm thick beryllium. A photomultiplier tube (PMT;

Received 4 September 2013

* Supported by NSFC (10978001, 11222327)

1) E-mail: jijianf@tsinghua.edu.cn

©2014 Chinese Physical Society and the Institute of High Energy Physics of the Chinese Academy of Sciences and the Institute of Modern Physics of the Chinese Academy of Sciences and IOP Publishing Ltd

HAMAMATSU R877) is coupled to the detector below the NaI(Tl). With simulations using GEANT4, the detection efficiency versus photon energy is shown in Fig. 1, where the solid curve corresponds to the efficiency of detection, i.e. the probability of any amount of energy deposited in LaBr₃, and the dashed curve is for the efficiency of full energy deposit in LaBr₃. At 662 keV, the detection efficiency is 22% while the full energy efficiency is only 5%.




Fig. 1. Simulated quantum efficiency of the phoswich detector as a function of energy. The jumps in the dashed curve correspond to the K shell binding energies of Br and La, respectively.

The PMT anode is connected to a preamplifier with a discharge timescale of roughly 10 ns, whose output is fed into a digitizer (CAEN V1742, 12-bit) sampling at a rate of 1 GS/s with a buffer size of 1024. Thus, every event recorded by the digitizer has a length of 1024 ns in time. The PMT was operated at a negative high voltage of 1080 V. All the experiments were conducted at room temperature.

3 System response and calibration

The light output of the scintillation crystal is expected to be an exponential curve,

$$I = I_0 e^{-\frac{t}{\tau}}, \quad (1)$$

where τ is the decay time of the scintillator and I_0 is the light yield per unit time at the very beginning. Such a signal will go through the PMT and the preamplifier and convolve with their responses. We define their overall response as the system response R_s . The system response can be measured by illuminating the PMT cathode with a pulsed light source whose pulse width is much smaller than the response time. A nano-second ultraviolet (UV) LED pulser (HORIBA NanoLED-265) is used to generate light pulses, with a peak wavelength of 268 nm and a full width at half maximum (FWHM) less than 1.2 ns.

The LED light is injected to the central position on the PMT cathode and the system response R_s is measured by the digitizer, see Fig. 2. The response shows a variability of less than 3.5% in FWHM with the change of PMT high voltage, illumination location and the light luminosity, which is negligible for the analysis.




Fig. 2. The measured system response R_s using a nano-second UV LED.

The output signal from the preamplifier for an X-ray event, S_o is therefore a convolution of the light decay with the system response,

$$S_o = IR_s. \quad (2)$$

4 Measurements and results

A radioactive ¹³⁷Cs source (source activity $\sim 10 \mu\text{Ci}$) is used to produce 662 keV gamma-rays. A measurement of FWHM versus amplitude for signals from the preamplifier is shown in Fig. 3. Each event is a point on




Fig. 3. FWHM versus amplitude (left) for measured events with a ¹³⁷Cs at 662 keV distributions. The right panel shows the distribution of FWHM and the two peaks are events with energy deposit in LaBr₃ and NaI, respectively. The dashed lines define two regions (80–90 ns and 260–320 ns, respectively) for events selection.

the left panel. The pulse width distribution is shown on the right, with two peaks corresponding to energy deposit in LaBr₃ and NaI, respectively (~ 85 ns vs. ~ 290 ns). There are some events with FWHM in between the LaBr₃ and NaI peaks, lying on a curve on the FWHM versus peak diagram, which are Compton events with energy deposit in both layers.

The energy spectra for pure LaBr₃ or NaI events are shown in Fig. 4. The LaBr₃ events are selected for those with FWHM in the range of 80–90 ns, and the NaI events are those in the 260–320 ns range. The channel-energy calibration is done using radioactive sources ⁶⁰Co(1173.2 keV, 1332.5 keV) and ²⁴¹Am(59.5 keV). The peak in Fig. 4 is the full energy deposit from the 662 keV gamma-rays and the broad continuum on its left is from Compton events with partial energy deposit.




Fig. 4. Energy spectra of ¹³⁷Cs for pure LaBr₃:Ce (top) or NaI(Tl) (bottom) events.

For events with FWHM between 90 ns and 260 ns, which have energy deposit in both layers, their signal waveforms are fitted using the following function

$$I = I_{10} e^{-\frac{t}{\tau_1}} + I_{20} e^{-\frac{t}{\tau_2}}, \quad (3)$$

where the subscription 1 and 2 denote LaBr₃ or NaI, respectively. The decay time is fixed at their theoretical values, $\tau_1 = 16$ ns and $\tau_2 = 250$ ns. There are only

two independent parameters in the fit, I_{10} and I_{20} , corresponding to the energy deposit in each crystal. A typical fitting result is shown in Fig. 5, where the two components are clearly disentangled. The total energy of the events can thus be reconstructed by adding up their energies in each crystal.

The reconstructed energy spectra for those “double deposit events” are shown in Fig. 6. The dashed curve is for pure LaBr₃ events and the solid curve is for the spectrum after adding up double deposit events. As one can see, the photopeak at 662 keV is higher but the energy resolution is slightly degraded due to worse resolution of NaI than LaBr₃. By including the double deposit events, the efficiency for full energy detection is increased by a factor of 1.9, consistent with the estimate using GEANT4 simulation.




Fig. 5. A typical fitting result of an event with energy deposit in both crystals. The measured and model (total, LaBr₃, and NaI) waveforms are shown on top and the data to model ratio is shown on bottom.




Fig. 6. Energy spectrum of ¹³⁷Cs consisting of pure LaBr₃ events (dashed) and by adding up double deposit events (solid).

5 Conclusion

In this work we demonstrate that Compton events in a LaBr_3/NaI phoswich detector with energy deposit in both $\text{LaBr}_3:\text{Ce}$ and $\text{NaI}(\text{Tl})$ can be recognized and the full energy of the incident photon can be reconstructed

using a digitizer. This improvement in readout can simply increase the efficiency of full energy detection by a factor of 2. This technique can be used in the energy range in which Compton events are significant, and may broaden the scientific return of phoswich detectors in X-ray and gamma-ray astronomy.

References

- 1 Frontera F, Costa E, Dal Fiume D, Feroci M, Nicastro L, Orlandini M, Palazzi E, Zavattini G. *Astronomy & Astrophysics Supplement Series*, 1997, **122**: 357–369
- 2 Shah K S, Glodo J, Klugerman M, Moses W W, Derenzo S E, Weber M J. *Nuclear Science, IEEE Transactions on*, 2003, **50**: 2410–2413
- 3 Mazumdar I, Kumar G A, Gothe D A, Manchanda R K. *Nucl. Instrum. Methods A*, 2010, **623**: 995–998
- 4 Manchanda R K. *Adv. Space Res.*, 2011, **47**: 30–36
- 5 LI H, JI J, FENG H, ZHANG Z, HAN D. *Nucl. Instrum. Methods A*, 2012, **671**: 24–28
- 6 Teegarden B J et al. *The Astrophysical Journal*, 2005, **621**: 296–300
- 7 Mallozzi R S, Pacias W S, Pendleton G N et al. *The Astrophysical Journal*, 1995, **454**: 597–603